

PRP-3, 12-20 µm Bulk Resin

Specification Sheet (Part/REF # 79702)



Pack your own column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Bulk Resin
Particle Size	12-20 µm
Packing Material Type	PSDVB
Pore Size	300 Å
Max Pressure	5,000 psi
Restoration	Run a gradient of 100% water to 100% acetonitrile (Repeat three times)
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79702>

Spec. Version B