

PRP-X500, 12-20 µm Bulk Resin

Specification Sheet (Part/REF # 79595)



Pack your own column.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X500
Format	Bulk Resin
Particle Size	12-20 µm
USP Classification	USP L23
Packing Material Type	Poly(methacryl amido-propyl trimethyl ammonium chloride)
Pore Size	Superficially Porous
Particle Surface Area (m²/cm³)	1.4
Pore Volume (cm³/g)	0.0063
Max Pressure	5,000 psi
Exchange Capacity	1.6 meq/gm
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79595>

Spec. Version B