

PRP-X600, 12-20 µm Bulk Resin

Specification Sheet (Part/REF # 79598)



Superficially porous, polymeric anion exchange coating for weak anion exchange (WAX) purification of nucleic acids such as single stranded/double stranded RNA and DNA, peptides, and proteins.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X600
Format	Bulk Resin
Particle Size	12-20 µm
Packing Material Type	PSDVB/Poly (dimethyl amido-propyl methacrylamide)
Pore Size	Superficially Porous
Particle Surface Area (m²/cm³)	1.4
Pore Volume (cm³/g)	0.0046
Exchange Capacity	1.6 meq/gm
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79598>

Spec. Version B