

22s gauge, Large Hub RN Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 7779-03)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

| | |
|----------------|--------------------|
| Gauge | 22s gauge |
| Point Style | 2: Curved, Beveled |
| Needle Length | 2.00 in (50.8 mm) |
| Autoclavable | Yes |
| Sales Quantity | 6 needles per pack |
| Hub | Large RN |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7779-03>

Spec. Version B