

25s gauge, Kel-F Hub Needle, 1.97 in, point style 3, 6/PK

Specification Sheet (Part/REF # 90149)



Syringes for U6K valve come with a 25s gauge needle specifically designed for this valve.

Product Specifications

Gauge	25s gauge
Point Style	3: Blunt
Needle Length	1.97 in (50.0 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Kel-F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/90149>

Spec. Version B