

# 18 gauge, Kel-F Hub Needle, 2.65 in, point style 3, 1/PK

Specification Sheet (Part/REF # 90432)



The VDRL Antigen Syringe is designed to simplify the dispensing steps for VDRL antigen testing. The test is a microscopic nontreponemal flocculation test to be used for the detection of reagin.

## Product Specifications

|                  |                   |
|------------------|-------------------|
| Gauge            | 18 gauge          |
| Point Style      | 3: Blunt          |
| Needle Length    | 2.65 in (67.3 mm) |
| Autoclavable     | Yes               |
| Gas Sterilizable | Yes               |
| Sales Quantity   | 1 needle per pack |
| Hub              | Kel-F             |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/90432>

Spec. Version B