

# 23 gauge, Large Hub RN Needle, custom length (0.5 to 12 in), point style 5 or AS, 6/PK

Specification Sheet (Part/REF # 7732-06)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	23 gauge
Point Style	5 or AS
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7732-06>

Spec. Version B