

26 gauge, Large Hub RN Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 7779-04)



Large Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 250 µL and larger, or SampleLock syringes. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	26 gauge
Point Style	2: Curved, Beveled
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7779-04>

Spec. Version B