

**PEEK 0.062 inch OD/ 0.030 inch ID,
Kel-F Hub Needle, 2 in, point style
3, 6/PK**

Specification Sheet (Part/REF # 8649-01)



Kel-F and Removable Needles using standard HPLC grade PEEK tubing.

Product Specifications

Gauge	0.062 inch OD/ 0.030 inch ID
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Kel-F
Needle Material	PEEK

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/8649-01>

Spec. Version B