

# 16 gauge, Kel-F Hub Needle, 2 in, point style 3, 6/PK

Specification Sheet (Part/REF # 90516)



Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available. These needles are available from 10 - 31 gauge.

## Product Specifications

Gauge	16 gauge
Point Style	3: Blunt
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Kel-F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/90516>

Spec. Version B