

22s gauge, Kel-F Hub Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 90138)



Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available. These needles are available from 10 - 31 gauge.

Product Specifications

| | |
|----------------|--------------------|
| Gauge | 22s gauge |
| Point Style | 2: Curved, Beveled |
| Needle Length | 2.00 in (50.8 mm) |
| Autoclavable | Yes |
| Sales Quantity | 6 needles per pack |
| Hub | Kel-F |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/90138>

Spec. Version B