

# 22s gauge, Metal Hub Chromatography Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 90051)



Metal (N) hub needles (nickel plated brass) are compatible with PTFE Luer Lock (TLL) syringes and plastic Luer Tip (LT) or Luer Lock (TLL) connectors. These 304 stainless steel needles are packaged in convenient 6 packs and are available from 10-30 gauge. Custom needle lengths and point styles are available.

## Product Specifications

|                |                    |
|----------------|--------------------|
| Gauge          | 22s gauge          |
| Point Style    | 2: Curved, Beveled |
| Needle Length  | 2.00 in (50.8 mm)  |
| Autoclavable   | Yes                |
| Sales Quantity | 6 needles per pack |
| Hub            | Metal (N)          |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/90051>

Spec. Version B