

25 gauge, Kel-F Hub Needle, custom length (0.375 to 12 in), point style 5 or AS, 6/PK

Specification Sheet (Part/REF # 7752-16)



Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available. These needles are available from 10 - 31 gauge.

Product Specifications

| | |
|------------------|--------------------|
| Gauge | 25 gauge |
| Point Style | 5 or AS |
| Needle Length | 0.75 in (19.0 mm) |
| Autoclavable | Yes |
| Gas Sterilizable | Yes |
| Sales Quantity | 6 needles per pack |
| Hub | Kel-F |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7752-16>

Spec. Version B