

25s gauge Tapered Point, Kel-F Hub Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 90152)



Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available. These needles are available from 10 - 31 gauge.

Product Specifications

Gauge	25s gauge
Point Style	2: Curved, Beveled
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Kel-F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/90152>

Spec. Version B