

# 26s gauge, Metal Hub Chromatography Needle, 2 in, point style 2, 6/PK

Specification Sheet (Part/REF # 90053)



Metal (N) hub needles (nickel plated brass) are compatible with PTFE Luer Lock (TLL) syringes and plastic Luer Tip (LT) or Luer Lock (TLL) connectors. These 304 stainless steel needles are packaged in convenient 6 packs and are available from 10-30 gauge. Custom needle lengths and point styles are available.

## Product Specifications

Gauge	26s gauge
Point Style	2: Curved, Beveled
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Metal (N)

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/90053>

Spec. Version B