

# 33 gauge, Metal Hub Needle, custom length (0.4 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7747-01)



Metal (N) hub needles (nickel plated brass) are compatible with PTFE Luer Lock (TLL) syringes and plastic Luer Tip (LT) or Luer Lock (TLL) connectors. These 304 stainless steel needles are packaged in convenient 6 packs and are available from 10-30 gauge. Custom needle lengths and point styles are available.

## Product Specifications

Gauge	33 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Metal (N)

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7747-01>

Spec. Version B