

22 gauge, Large Hub RN Needle,
custom length (1.00 to 6.5 in), point
style 3T, 6/PK

Specification Sheet (Part/REF # 8646-01)



The PTFE coated needles are used for thin layer chromatography (TLC) and are coated with PTFE 3/4" up the needle from the tip.

Product Specifications

Gauge	22 gauge
Point Style	3T: Blunt PTFE Coated
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Large RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/8646-01>

Spec. Version B