

18 gauge, Small Hub RN Needle, custom length (0.5 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7804-06)

Small Gauge Removable Needle (RN) Hub Needle

This hub is used for needles between 27-34 gauge. Needle length is specified from the needle point to the start of the support sleeve.



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	18 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7804-06>

Spec. Version B