

# 21 gauge, Small Hub RN Needle, custom length (0.5 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7804-12)



## Needle Point Style 2

Sharp, beveled, curved, non-coring needle point recommended for piercing septa. Available gauges: 34-10.

## Needle Point Style 3

Blunt needle point for use with HPLC injection valves and for sample pipetting. Available gauges: 34-10.

## Needle Point Style 3T

Blunt needle point for thin layer chromatography (TLC) applications. Needles are coated with PTFE 3/4" up from the needle tip. Available gauges: 26s, 26, 22s, 22.

Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	21 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7804-12>

Spec. Version B