

# 22s gauge, Small Hub RN Needle, 2 in, point style 2, 6/Pack

Specification Sheet (Part/REF # 7758-03)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	22s gauge
Point Style	2: Curved, Beveled
Needle Length	2.00 in (50.8 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7758-03>

Spec. Version B