

22s gauge, Small Hub RN Needle,
custom length (1.00 to 6.5 in), point
style 3T, 6/PK

Specification Sheet (Part/REF # 8646-03)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

| | |
|----------------|-----------------------|
| Gauge | 22s gauge |
| Point Style | 3T: Blunt PTFE Coated |
| Needle Length | Custom Needle Length |
| Autoclavable | Yes |
| Sales Quantity | 6 needles per pack |
| Hub | Small RN |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/8646-03>

Spec. Version B