

# 23s-26s gauge, Small Hub RN Needle, 1.71 in, point style AS, 6/PK

Specification Sheet (Part/REF # 7785-01)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	23s-26s gauge
Point Style	AS: Conical, Autosampler
Needle Length	1.69 in (42.9 mm)
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7785-01>

Spec. Version B