


# 30 gauge, Small Hub RN Needle, custom length (0.5 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7803-07)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

Gauge	30 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/7803-07>

Spec. Version B