

33 gauge, Small Hub RN Needle, 1.5 in, point style 3, 6/Pack

Specification Sheet (Part/REF # 7762-06)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	33 gauge
Point Style	3: Blunt
Needle Length	1.50 in (38.1 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7762-06>

Spec. Version B