

**34 gauge, Small Hub RN Needle,
custom length (0.375, 0.5, 1.0, and
1.5 in), point style 3 or 4, 6/PK**

Specification Sheet (Part/REF # 207434)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	34 gauge
Point Style	3 or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/207434>

Spec. Version B