

25s gauge, Hubless Needle, custom length (0.4 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 22058-01)



Some custom applications require a cannula only without the needle hub attached.

Product Specifications

| | |
|----------------|----------------------|
| Gauge | 25s gauge |
| Point Style | 2, 3, or 4 |
| Needle Length | Custom Needle Length |
| Autoclavable | Yes |
| Sales Quantity | 6 needles per pack |
| Hub | No Hub |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/22058-01>

Spec. Version B