

# 26s gauge, Raw Tubing, Sold Per Inch

Specification Sheet (Part/REF # 21039A)



Hamilton offers raw cut stainless steel tubing by the inch.

## Product Specifications

|                |                      |
|----------------|----------------------|
| Gauge          | 26s gauge            |
| Needle Length  | Custom Needle Length |
| Autoclavable   | Yes                  |
| Sales Quantity | Sold Per Inch        |
| Hub            | No Hub               |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/21039a>

Spec. Version B