

33 gauge, Raw Tubing, Sold Per Inch

Specification Sheet (Part/REF # 21033A)



Hamilton offers raw cut stainless steel tubing by the inch.

Product Specifications

| | |
|------------------|----------------------|
| Gauge | 33 gauge |
| Needle Length | Custom Needle Length |
| Autoclavable | Yes |
| Gas Sterilizable | Yes |
| Sales Quantity | Sold Per Inch |
| Hub | No Hub |

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/21033a>

Spec. Version B