

**14 gauge, Hubless Needle, custom length (0.4 to 12 in), point style 2, 3, or 4, 6/PK**

Specification Sheet (Part/REF # 22014-01)



Some custom applications require a cannula only without the needle hub attached.

**Product Specifications**

Gauge	14 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	No Hub

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/22014-01>

**Spec. Version B**