

23 gauge, Hubless Needle, custom length (0.4 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 22023-01)



Some custom applications require a cannula only without the needle hub attached.

Product Specifications

Gauge	23 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	No Hub

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/22023-01>

Spec. Version B