

28 gauge, Hubless Needle, custom length (0.4 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 22028-01)



Some custom applications require a cannula only without the needle hub attached.

Product Specifications

Gauge	28 gauge
Point Style	2, 3, or 4
Needle Length	Custom Needle Length
Autoclavable	Yes
Sales Quantity	6 needles per pack
Hub	No Hub

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/22028-01>

Spec. Version B