

# Adjustable Volume SoftGrip Pipette

## 100-1000 µL

Specification Sheet (Part/REF # 55019-40)



Available in 0.2µL - 1 mL, SoftGrip’s isolated volume adjustment increases accuracy by removing parallax errors and preventing accidental volume adjustment, thanks to our next-generation dual indicator. Hamilton pipettors boast a slim design that allows for reaching the bottom of narrow tubes without the need for extended-length tips.

### Product Specifications

Volume Range	100 to 1000 µL
Channels	1
Color	Sky Blue
Accuracy at 100 Percent Volume	0.30 ±%
Accuracy at 10% Volume	1.60 ±%
Precision at 100 Percent Volume	0.12 ±%
Precision at 10% Volume	0.50 ±%
Increments	2.0 µL
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/pipettes/55019-40>

Spec. Version B