

# Adjustable Volume SoftGrip Pipette

## 0.2-2 $\mu$ L

Specification Sheet (Part/REF # 55019-30)



Available in 0.2 $\mu$ L - 1 mL, SoftGrip's isolated volume adjustment increases accuracy by removing parallax errors and preventing accidental volume adjustment, thanks to our next-generation dual indicator. Hamilton pipettors boast a slim design that allows for reaching the bottom of narrow tubes without the need for extended-length tips.

### Product Specifications

Volume Range	0.2 to 2 $\mu$ L
Channels	1
Color	Aqua
Accuracy at 100 Percent Volume	1.20 $\pm$ %
Accuracy at 10% Volume	8.00 $\pm$ %
Precision at 100 Percent Volume	0.60 $\pm$ %
Precision at 10% Volume	4.00 $\pm$ %
Increments	0.002 $\mu$ L
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/pipettes/55019-30>

Spec. Version B