

Fixed Volume SoftGrip Pipette 10 µL

Specification Sheet (Part/REF # 55019-03)



SoftGrip fixed volume pipettes come with a slim ejection sleeve made of chemically resistant, high-strength polymer, as well as a universal calibration key to reduce calibration time. Hamilton fixed volume pipettes are available in 5 µL - 1 mL.

Product Specifications

Volume Range	10 µL Fixed
Channels	1
Color	Purple
Accuracy at 100 Percent Volume	0.80 ±%
Accuracy at 10% Volume	N/A
Precision at 100 Percent Volume	0.40 ±%
Precision at 10% Volume	N/A
Increments	N/A
Autoclavable	Yes
Field Calibration	Yes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/pipettes/55019-03>

Spec. Version B