

# Fixed Volume SoftGrip Pipette 300 $\mu$ L

Specification Sheet (Part/REF # 55019-15)



SoftGrip fixed volume pipettes come with a slim ejection sleeve made of chemically resistant, high-strength polymer, as well as a universal calibration key to reduce calibration time. Hamilton fixed volume pipettes are available in 5  $\mu$ L - 1 mL.

## Product Specifications

|                                 |                   |
|---------------------------------|-------------------|
| Volume Range                    | 300 $\mu$ L Fixed |
| Channels                        | 1                 |
| Color                           | Brick Red         |
| Accuracy at 100 Percent Volume  | 0.40 $\pm$ %      |
| Accuracy at 10% Volume          | N/A               |
| Precision at 100 Percent Volume | 0.18 $\pm$ %      |
| Precision at 10% Volume         | N/A               |
| Increments                      | N/A               |
| Autoclavable                    | Yes               |
| Field Calibration               | Yes               |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/pipettes/55019-15>

Spec. Version B