

10 µL, Bulk, Low Retention, 1000 tips/bag

Specification Sheet (Part/REF # 68007-20)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

Product Specifications

Disposable Tip Volume	10 µL
Tip Length	Standard
Packaging	Bulk
Quantity Per Each	1000
Filtered	No
Low Retention	Yes
Sterile	No
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recomended Pipette Volume Ranges	2 µL, 5 µL, and 10 µL Pipettes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/pipettes/68007-20>

Spec. Version B