

# 200 µL, Racked, 96 tips/rack

Specification Sheet (Part/REF # 68008-11)



Hamilton pipette tip families deliver value, choice, and performance across a variety of formats and packaging options. Formats include standard, non-filter, extended length, hydrophobic filter barrier, non-sterile, sterilized, and macro volume.

## Product Specifications

Disposable Tip Volume	200 µL
Tip Length	Standard
Packaging	Racked
Quantity Per Each	960
Filtered	No
Low Retention	No
Sterile	No
Purity Level	Free of RNase/DNase, ATP, endotoxin, cytotoxic effects, and trace heavy metals
Recomended Pipette Volume Ranges	20 µL, 25 µL, 50 µL, 100 µL, and 200 µL Pipettes

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/pipettes/68008-11>

Spec. Version B