

# DuraCal pH buffer 4.01/7.00/10.01

## 500 ml

Specification Sheet (Part/REF # 238924)



Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

### Product Specifications

pH Accuracy	± 0.01 / ± 0.02 pH
Package	500 mL, each
Stability Months (after date of manufacturing)	± 0.01 pH: 24 / ± 0.02 pH: 60
Certificate	DAkkS
Storage Conditions	5 - 35 °C / 41 - 95 °F
pH Value	4.01/7.00/10.01

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/sensors/238924>

Spec. Version C