

# ORP buffer +475 mV, 500 ml

Specification Sheet (Part/REF # 238227)



Hamilton ORP buffers are stable for two years. You can always rely on them even if they have been used.

## Product Specifications

Package	500 mL
Redox Value mV	475 mV
mV Accuracy	±5 mV
Stability Months (after date of manufacturing)	24
Certificate	No
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/sensors/238227>

Spec. Version C