

ORP buffer +271 mV, 500 ml

Specification Sheet (Part/REF # 238228)



Hamilton ORP buffers are stable for two years. You can always rely on them even if they have been used.

Product Specifications

Package	500 mL
Redox Value mV	271 mV
mV Accuracy	±5 mV
Stability Months (after date of manufacturing)	24
Certificate	No
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/sensors/238228>

Spec. Version c