

ORP buffer +475 mV, 250 ml

Specification Sheet (Part/REF # 238322)



Hamilton ORP buffers are stable for two years. You can always rely on them even if they have been used.

Product Specifications

Package	250 mL
Redox Value mV	475 mV
mV Accuracy	±5 mV
Stability Months (after date of manufacturing)	24
Certificate	No
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/sensors/238322>

Spec. Version C