

# Conductivity standard 100 mS/cm

Specification Sheet (Part/REF # 238935)



Stability from 1 to 3 years respectively and conductivity ranges from 1.3  $\mu$ S/cm to 100 mS/cm

## Product Specifications

Accuracy at 25 °C	± 1 %
Value at 25 °C	100 mS/cm
Package	250mL glass bottle
Stability Months (after date of manufacturing)	36
Certificate	Hamilton
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/sensors/238935>

Spec. Version D