

Conductivity standard 1.3 μS/cm

Specification Sheet (Part/REF # 238973)



Stability from 1 to 3 years respectively and conductivity ranges from 1.3 μS/cm to 100 mS/cm

Product Specifications

Accuracy at 25 °C	± 1 %
Value at 25 °C	1.3 μS/cm
Package	250mL glass bottle
Stability Months (after date of manufacturing)	12
Certificate	DFM
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/sensors/238973>

Spec. Version C