

# Conductivity standard 1413 µS/cm, Basic Line

Specification Sheet (Part/REF # 238986)



Stability from 1 to 3 years respectively and conductivity ranges from 1.3 µS/cm to 100 mS/cm

## Product Specifications

Accuracy at 25 °C	± 1 %
Value at 25 °C	1,413 µS/cm
Package	500 mL Calpack bottle
Stability Months (after date of manufacturing)	18
Certificate	DFM
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/sensors/238986>

Spec. Version C