

Conductivity standard 147 µS/cm, Basic Line

Specification Sheet (Part/REF # 238985)



Stability from 1 to 3 years respectively and conductivity ranges from 1.3 µS/cm to 100 mS/cm

Product Specifications

Accuracy at 25 °C	± 1 %
Value at 25 °C	147 µS/cm
Package	500 mL Calpack bottle
Stability Months (after date of manufacturing)	18
Certificate	DFM
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/sensors/238985>

Spec. Version C