

DuraCal pH buffer 3.06, 500 ml

Specification Sheet (Part/REF # 238274)



Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

Product Specifications

pH Accuracy	± 0.02 pH
Package	500 mL
Stability Months (after date of manufacturing)	60
Certificate	Hamilton
Storage Conditions	5 - 35 °C / 41 - 95 °F
pH Value	3.06

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/sensors/238274>

Spec. Version C