

# DuraCal pH buffer pH 10.01, 3 x 500 ml

Specification Sheet (Part/REF # 238923)



Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

## Product Specifications

pH Accuracy	± 0.02 pH
Package	3 x 500 mL
Stability Months (after date of manufacturing)	60
Certificate	DAkkS
Storage Conditions	5 - 35 °C / 41 - 95 °F
pH Value	10.01

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/sensors/238923>

Spec. Version C