

PB600-1 Dispenser

Specification Sheet (Part/REF # 83700)



The VDRL Antigen Syringe is designed to simplify the dispensing steps for VDRL antigen testing. The test is a microscopic nontreponemal flocculation test to be used for the detection of reagin.

Product Specifications

Compatible with 700 and 1700 series syringes in the following volumes: 10 μ L, 25 μ L, 50 μ L, 100 μ L, 250 μ L, 500 μ L, 1 mL, and 2.5 mL

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/other-laboratory/83700>

Spec. Version B