

# Micro Pipette Disposable Tip, PTFE

## 100 pack

Specification Sheet (Part/REF # 84255)



Microsyringe Pipettes have a special 3/4" electrotapered needle for easy attachment of PTFE tips. Your sample only wets the disposable PTFE tips, eliminating sample cross-contamination. The spring-loaded plunger and adjustable stops ensure excellent reproducibility from 0.2 to 3  $\mu$ L.

### Product Specifications

Sales Quantity	100 tips per pack
Specifications are subject to change without notice <a href="https://www.hamiltoncompany.com/syringes/84255">https://www.hamiltoncompany.com/syringes/84255</a>	
Spec. Version B	