

5 mL, Model 1005 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 9882-45)



Product Specifications

Volume	5 mL
Syringe Type	Longlife Gastight
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Plunger Style	Side Post Mounting (XL)
Fluid Path	Borosilicate glass, PTFE, UHMWPE
Autoclavable	No
Gas Sterilizable	Yes
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.406 inches (10.3 mm)
Barrel Outer Diameter	0.530 inches (13.5 mm)
Maximum Pressure	200 psig (13.8 bar)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/9882-45>

Spec. Version B